

Candidate's name:.....

Signature:.....COMBINATION:.....

P525/2  
CHEMISTRY  
Paper 2  
(Practical)  
July/Aug.2025  
2 hours

THE CHEMISTRY DEPARTMENT

CHEMISTRY

S5

Paper 2  
(Practical)

2 hours

INSTRUCTIONS TO CANDIDATES:

*This paper consists of **one compulsory** examination item.*

*Answers to this item are to be written in the spaces provided in this paper. Use **blue** or **black** ink. Any work done in pencil **except** drawings and graph will **not** be scored.*

***All** working **must** be clearly shown. Graph paper will be provided.*

*Mathematical tables and silent non-programmable scientific calculators may be used.*

*You are **not** allowed to use reference books.*

*Candidates are advised to carefully read the item, make sure they have all the apparatus and chemicals they may need and then **plan** appropriately before starting.*

TURNOVER

## ITEM


A school chemistry lab uses a sodium hydroxide solution of concentration at approximately 40 g/dm<sup>3</sup> for acid-qualitative analysis. Over time, the solution absorbs carbon dioxide from air, forming sodium carbonate and reducing its concentration. The lab technician suspects the current solution may no longer meet the required concentration.

Accurate concentration is critical for an upcoming test on unknown carboxylic acids. An incorrect sodium hydroxide concentration could lead to faulty interpretations.

To verify its suitability, there is need to standardise the sodium hydroxide solution using a primary standard solution of oxalic acid (H<sub>2</sub>C<sub>2</sub>O<sub>4</sub>·2H<sub>2</sub>O), prepared by dissolving 6.90–7.10 g in 250 cm<sup>3</sup> of distilled water.

You will then determine whether the sodium hydroxide solution's concentration is approximately 40g per litre and advise whether it can be used as-is, or must be discarded. The outcome will ensure the lab's acid-base experiments and the students' results remain accurate and reliable.

Sodium hydroxide reacts with oxalic acid according to the following equation;



You are provided with the following;

- 7.5 g of Oxalic acid (H<sub>2</sub>C<sub>2</sub>O<sub>4</sub>·2H<sub>2</sub>O)
- distilled water, sufficient for solution preparation and rinsing
- 100cm<sup>3</sup> of GA1 which is the solution of sodium hydroxide used by the laboratory.
- phenolphthalein indicator
- 250.0 cm<sup>3</sup> volumetric flask
- Other laboratory apparatus suitable for a titration experiment.

## Task;

As a chemistry learner;

- [illegible]

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

(b) Carry out the experiment and record your results.

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

(c) Interpret and analyse your results.

.....

.....

.....

.....

This image shows a full page of primary-ruled paper. It features multiple horizontal rows, each defined by two parallel dotted lines. The rows are evenly spaced across the entire page, providing a guide for handwriting practice. There are no margins, text, or other markings present.

(d) What conclusion and recommendation can you give basing on your results?

.....

.....

.....

**END**